CLINICAL EVALUATION REPORT

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Product: COOLING GEL SHEET

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Document No.: QD/QB-730-01-11

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Version: A0

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Date: 2024-03-06

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Author: He Jian Qiang

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Reviewed by: Rong Cheng

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Approved by: Chen Ning

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CLINICAL EVALUATION REPORT Cooling Gel Sheet 冷贴临床评估报告	Document No.: QD/QB-730-01-11 文件编号: QD/QB-730-01-11
	Version: A/0 版本: A/0
	Date: 2024-03-06 日期: 2024-03-06

Table of Content

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	Page 页码
1. General details 总述	3
2. Description of the device and its intended application 器械描述和预期用途	3
3. Intended therapeutic and/or diagnostic indications and claims 预期治疗和/或诊断说明和要求	4
4. Context of the evaluation and choice of clinical data types 评估背景和临床数据类型的选择	4
5. Summary of the clinical data and appraisal 总结临床数据和评价	5
6. Data analysis 数据分析	7
6.1. Performance 性能	7
6.2. Safety 安全	8
6.3. Product Literature and Instructions for Use 产品文献和使用说明	9
7. Conclusions 结论	10

CLINICAL EVALUATION REPORT Cooling Gel Sheet 冷贴临床评估报告	Document No.: QD/QB-730-01-11 文件编号: QD/QB-730-01-11
	Version: A/0 版本: A/0
	Date: 2024-03-06 日期: 2024-03-06

1. General details 总述

State the proprietary name of the device and any code names assigned during device development.

Identify the manufacturer(s) of the device.

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Name of the device: Cooling Gel Sheet

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Manufacturer: Foshan Biours Biosciences Co., Ltd.

Description of the device and its intended application 器械描述和预期用途

Provide a concise physical description of the device, cross-referencing to relevant sections of the manufacturer's technical information as appropriate. The description should cover information such as: materials, including whether it incorporates a medicinal substance (already on the market or new), tissues, or blood products; the device components, including software and accessories; mechanical characteristics; and others, such as sterile vs. non-sterile, radioactivity etc.

State the intended application of the device, single use/reusable; invasive/noninvasive; implantable; duration of use or contact with the body; organs, tissues or body fluids contacted by the device.

Describe how the device achieves its intended purpose.

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Composed of kaolin, sodium polyacrylate, aluminum dihydroxyaminoacetate, disodium ethylenediaminetetraacetate, blue/green/red/orange pigments, glycerin, purified water, jacket and fixture. Through the evaporation of water, it absorbs heat and takes away the heat generated when the human body has a fever, thereby achieving the effect of cooling and relieving fever. Does not contain ingredients that exert pharmacological, immunological or metabolic effects, and does not contain ingredients listed in the appendix "List of Prohibited Ingredients for Certain Class I Medical Device Products". Non-sterile products.

CLINICAL EVALUATION REPORT Cooling Gel Sheet 冷贴临床评估报告	Document No.: QD/QB-730-01-11 文件编号: QD/QB-730-01-11
	Version: A/0 版本: A/0
	Date: 2024-03-06 日期: 2024-03-06

2. Intended therapeutic and/or diagnostic indications and claims 预期治疗和/或诊断说明和要求

State the medical conditions to be treated, including target treatment group and diseases. Outline any specific safety or performance claims made for the device.

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For local cooling of patients with fever. For use on intact skin only. There are no special safety or performance requirements for this device.

3. Context of the evaluation and choice of clinical data types 评估背景和临床数据类型的选择

Outline the developmental context for the device. The information should include whether the device is based on a new technology, a new clinical application of an existing technology, or the result of incremental change of an existing technology. The amount of information will differ according to the history of the technology. Where a completely new technology has been developed, this section would need to give an overview of the developmental process and the points in the development cycle at which clinical data have been generated. For long-standing technology, a shorter description of the history of the technology (with appropriate references) could be used. Clearly state if the clinical data used in the evaluation are for an equivalent device. Identify the equivalent device(s) and provide a justification of the equivalency, cross-referenced to the relevant non-clinical documentation that supports the claim.

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This device technology originated from Japanese cataplast technology and has a long history of application in Japan. For conventional technology. The clinical data used in the evaluation of this device were obtained from equivalent devices.

State the Essential Requirements relevant to the device in question, in particular, any special design features that pose special performance or safety concerns (e.g. presence of medicinal, human or animal components) that were identified in the device risk management documentation and that required assessment from a clinical perspective.

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医疗器械的基本要求是安全有效。The basic requirement for medical devices is safety and effectiveness

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Cooling Gel Sheets have undergone necessary safety assessments from a clinical perspective. They are directly attached to human skin and have been identified in risk management documents. Cooling Gel Sheet products have been inspected in accordance with design requirements and are all qualified. (See test report)

Outline how these considerations were used to choose the types of clinical data used for the evaluation. Where published scientific literature has been used, provide a brief outline of the searching/retrieval process, cross-referenced to the literature search protocol and reports.

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Factors that need to be considered depending on the characteristics of the cooling gel sheets include the age and physical condition of the user. Age factor: Infants and children should be used under adult supervision. Physical condition factors: Do not apply it to wounds, mucous membranes, eczema and other skin injuries. Clinical data that exclude the above factors should be excluded from clinical evaluation.

4. Summary of the clinical data and appraisal 总结临床数据和评价

Provide a tabulation of the clinical data used in the evaluation, categorized according to whether the data address the performance or the safety of the device in question. (Note: many individual data sets will address both safety and performance.) Within each category, order the data according to the importance of their contribution to establishing the safety and performance of the device and in relation to any specific claims about performance or safety. Additionally, provide a brief outline of the data appraisal methods used in the evaluation, including any weighting criteria, and a summary of the key results.

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After statistical analysis for evaluation data, all data are in scope of standard requirements.

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Chart of the cooling patches clinical evaluation data

Include full citations for literature-based data and the titles and investigation codes (if relevant) of any clinical investigation reports.

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Cross-reference the entry for each piece of data to its location in the manufacturer's technical documentation.

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见生产记录和检验记录。See record of production and inspection record.

5. Data Analysis 数据分析

5.1. Performance 性能

Provide a description of the analysis used to assess performance.

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Cooling performance: Take the sample together with the intact packaging and place it at $23^{\circ}\text{C}\pm 1^{\circ}\text{C}$ and relative humidity of 40% to 70% for more than 8 hours for pretreatment. In a constant temperature and humidity box with a temperature maintained at $40^{\circ}\text{C} \pm 1^{\circ}\text{C}$ and a relative humidity of $50\% \pm 10\%$, four thermocouple temperature probes are installed. One probe is used to measure the real-time temperature of the space in the constant temperature and humidity box (hereinafter referred to as Thermocouple A), and the other three thermocouple probes are wrapped in the product and used to measure the real-time temperature of the sample under simulated use (hereinafter referred to as

CLINICAL EVALUATION REPORT Cooling Gel Sheet 冷贴临床评估报告	Document No.: QD/QB-730-01-11 文件编号: QD/QB-730-01-11
	Version: A/0 版本: A/0
	Date: 2024-03-06 日期: 2024-03-06

Thermocouple B, Thermocouple C, Thermocouple D). When thermocouple A measures the real-time temperature of the space in the constant temperature and humidity box to 40 °C, take the pretreated products and remove the anti-adhesive layer, and make the gel surface tightly wrap around the thermocouple in the constant temperature and humidity box. B. Thermocouple C, Thermocouple D, start timing, and read the temperature measured by thermocouple A at 10min, 20min, 30min, 40min, 50min, 60min, 90min, 120min, 150min, 180min, 210min, and 240min respectively. is Tai, the temperatures measured by thermocouple B, thermocouple C, and thermocouple D are Tbi, Tci, and Tdi respectively.

Identify the datasets that are considered to be the most important in contributing to the demonstration of the overall performance of the device and, where useful, particular performance characteristics. Outline why they are considered to be pivotal and how they demonstrate the performance of the device collectively (e.g. consistency of results, statistical significance, clinical significance of effects).

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Data on the overall performance and specific characteristics of this device are cooling properties. The main function of cooling gel sheet is to cool down.

5.2. Safety 安全

Describe the total experience with the device, including numbers and characteristics of patients exposed to the device; and duration of follow-up of device recipients.

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Tear open the packaging bag along the gap, take out the patch, uncover the protective film, and apply the gel surface directly to the skin area that needs cooling and cleansing. It can be cut appropriately as needed before use.

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1. Do not apply it to wounds, mucous membranes, eczema and other skin injuries.

2. Disposable external products, please do not reuse them.

3. Please clean the application area and keep it dry before use to prevent sweat, hair, etc. from affecting the viscosity of the product.

4. If local skin discomfort occurs during use, stop using it immediately. In severe cases, please consult a physician.

5. Infants and children should use it under adult supervision to prevent misuse and ingestion.

6. This product is used to assist cooling. If the fever persists, you should receive medical treatment.

Provide a summary of device-related adverse events, paying particular attention to serious adverse events.

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No device-related adverse events get involved.

Provide specific comment on whether the safety characteristics and intended purpose of the device requires training of the end-user.

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The device requires end-users to operate by instruction book, no need training.

5.3. Product Literature and Instructions for Use 产品文献和使用说明

State whether the manufacturer's proposed product literature and Instructions for Use are consistent with the clinical data and cover all the hazards and other clinically relevant information that may impact on the use of the device.

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Statement: the product literature and instructions for use of the product mentioned above are consistent with the clinical data and cover all the hazards and other clinically relevant information that may impact on the use of the device.

CLINICAL EVALUATION REPORT Cooling Gel Sheet 冷贴临床评估报告	Document No.: QD/QB-730-01-11 文件编号: QD/QB-730-01-11
	Version: A/0 版本: A/0
	Date: 2024-03-06 日期: 2024-03-06

6. Conclusions 结论

Outline clearly the conclusions reached about the safety and performance of the device from the evaluation, with respect to the intended use of the device. State whether the risks identified in the risk management documentation have been addressed by the clinical data.

For each proposed clinical indication state whether:

the clinical evidence demonstrates conformity with relevant Essential Requirements;

the performance and safety of the device as claimed have been established;

and the risks associated with the use of the device are acceptable when weighed against the benefits to the patient.

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The evaluation of the clinical data involves risks identified in risk management documentation. The conclusion is when the end-user is operating the device, safety and risks associated with the use of the device are acceptable, with respect to the intended use of the device.